

# General Purpose Infusion Systems

Description	Product Code	REF Number	Qty.
-------------	--------------	------------	------

## Outlook® 100 Safety Infusion System



The Outlook 100 Safety Infusion System is designed to help prevent IV medication errors. The Outlook 100 Infusion Device uses DoseGuard™ technology to establish and assign dose limits to critical medications. The DoseGuard software notifies the clinician when these hospital-defined limits are exceeded during manual programming and titration.

Outlook 100 Infusion Device	620-100	620-100	1
Syringe Holder	FZ0054	FZ0054	24
WinDock™ Programming Kit	FZ501630	FZ501630	1

## Outlook 200 Safety Infusion System



The Outlook 200 Safety Infusion System incorporates both the DoseGuard dose limit technology found in the Outlook 100 as well as the patented DoseScan™ technology. The DoseScan system uses bar codes to help ensure the *Right Patient* receives the *Right Medication* at the *Right Dose* by the authorized clinician at the point of care. The DoseScan feature helps prevent IV medication errors by providing patient matching and automated programming capabilities.

Outlook 200 Infusion Device	620-200	620-200	1
Syringe Holder	FZ0054	FZ0054	24
WinDock™ Programming Kit	FZ501630	FZ501630	1
Hand-Held Bar Code Scanner	FZ0098	FZ0098	1

# General Purpose Infusion Systems

Description	Product Code	REF Number	Qty.
-------------	--------------	------------	------

## Outlook® 300 Safety Infusion System



The Outlook 300 Safety Infusion System incorporates all the features and technology found in the Outlook 200 as well as the patented DoseCom™ feature. The DoseCom feature in the Outlook 300 allows the transmission of real-time infusion data across the hospital's existing 802.11a, b, or g wireless network. The Outlook 300 is also compatible with the DoseTrac™ Infusion Management Software for evaluating the safety, efficiency and efficacy of your IV therapy.

Outlook 300 Infusion Device	620-300	620-300	1
Syringe Holder	FZ0054	FZ0054	24
WinDock™ Programming Kit	FZ501630	FZ501630	1

## Outlook 400 Safety Infusion System



The Outlook 400 Safety Infusion System incorporates all the features and technology found in the Outlook 300 plus the ability to utilize bi-directional wireless data to provide a method to customize communications for your institutions infusion therapy needs. The Outlook 400 is also compatible with the DoseTrac Infusion Management Software for evaluating the safety, efficiency and efficacy of your IV therapy.

Outlook 400 Infusion Device	620-400	620-400	1
Syringe Holder	FZ0054	FZ0054	24
WinDock Programming Kit	FZ501630	FZ501630	1

## DoseTrac Infusion Management Software



The DoseTrac Infusion Management Software (IMS) provides the ability to capture transmitted data from the Outlook 300 and Outlook 400 infusion devices to help you proactively evaluate your IV therapy program. All members of the infusion team can benefit from the real-time views and custom reporting capabilities of the DoseTrac software. Optional scheduled reporting capabilities allow authorized users to define custom report formats for automatic delivery.

DoseTrac Infusion Management S/W	FZ620DTIMS	FZ620DTIMS	1
DoseTrac IMS w/ Report Scheduler	FZ620DTWRS	FZ620DTWRS	1